

EXOGUARD

The Ultimate Team Sports Protection



Order Form & Video Instructions

Access the Order Form here >



https://www.thuasneusa.com/exoguard_fill-out/

Watch Scanning Video >



<https://bit.ly/2FiHgY3>

Watch Casting Video >



<https://bit.ly/2Fjfqeb>

For additional information, speak to your local Townsend/Thuasne USA representative or call our fabrication center to speak with a specialist about proper casting and scanning protocols.



EXOGUARD



5545 St. Jacques O. Montreal, QC H4A 2E3

O: 514. 369. 3311

Toll Free / Sans Frais : 1.877.222.3311 🇨🇦 Fax Toll Free / Telec Sans Frais : 1.877.527.3311

www.bmiortho.com




Ultimate Protection

BMI
BERRETTA MEDICAL INC.

EXOGUARD

Ultimate Protection

Combining the most advanced composite technology with Townsend's tradition of perfect fit and function, EXOGUARD provides the ultimate protection for athletes.



Ultimate Precision

- The Townsend Motion joint mirrors knee kinematics providing continuous brace-to-leg contact, enhancing functional control.
- The unique tibial shell design and vacuum lamination process provides superior rotary control.
- No migration and no pistoning thanks to the combination of a biomechanically correct hinge and our patented suspension strap technology.

Ultimate Strength

- Advanced composite technology creates an extremely durable yet lightweight structure.
- Independent lab results prove EXOGUARD has superior durability in lateral impact testing.
- HD lamination option specially designed for high impact sports.
- Available with aircraft grade aluminum or stainless steel hinges.



Ultimate Fit

- EXOGUARD is a true custom knee brace, molded to a positive model of the leg, obtained either from a cast or a scan.
- Intimate shell to skin surface contact provides superior control and protection.